



Touch Screen Electronic Analytical Balance

Product Description

Feature:

- Touch Screen Electronic Analytical Balance
- 4.3-inch color touch screen
- Powerful functions like weighing, parts counting, percent weighing, check weighing, animal weighing, density weighing, etc.
- Electromagnetic balance sensor
- External/ Internal calibration
- Supports PC connection. Wireless Bluetooth printing is optional
- The sealed arc-shaped cover is made of acrylic material, beautiful appearance and easy to use



DIMENSION





Cat#	Weighing Range	Calibration	Readability	Stabilization Time	Pan Size	Dimension / N.W.
01-BD104EX	0-100g	External	0.1mg	≤4S	Φ85mm	350*215*340mm/8KG
01-BD204EX	0-200g					
01-BD224EX	0-220g					
01-BD304EX	0-300g					
01-BD104NX	0-100g	Internal				
01-BD204NX	0-200g					
01-BD224NX	0-220g					
01-BD304NX	0-300g					
01-BD103EX	0-100g	External	1mg			
01-BD203EX	0-200g					
01-BD263EX	0-260g					
01-BD303EX	0-300g					
01-BD503EX	0-500g					
01-BD103NX	0-100g	Internal				
01-BD203NX	0-200g					
01-BD263N1	0-260g					
01-BD303NX	0-300g					
01-BD503NX	0-500g					
Standard Accessories	Calibration Weight (for external calibration), Adapter, Manual, Dust Cover					