

# Major Science BluView Transilluminator MBE-200

## Product Introduction

The BluView Transilluminator is packed with powerful 470nm blue LED light to help you enhance the gel images while maintain environmental safe to you and the mother Earth. With the gel imaging size 15.4x15.4cm, you could fit any midi to small size gel on top of the transilluminator. Attached a 580nm filter lid allows you to view the sample within any angle available, adjustable hinge prevents any glare or reflection might affect your visibility.

The BluView Transilluminator is also compatible with our imaging systems such as Digimage system and SmartView Simple Imager systems.

The BluView Transilluminator is also compatible with our imaging systems such as Digimage system and SmartView Simple Imager systems.

## Features

- Light Weight, easy to carry in between labs
- Aluminium alloy casing design
- Low heat dispersion
- Energy saving product
- 470nm blue light for safe human contact



## Flexible

580nm filter lid for viewing at any angle available without glare or reflection.

## Optimized

470nm LED Blue light, providing best performance with most safe dye sensitivity: 0.46 ng

