

Midi plus-1 Horizontal Electrophoresis System, ME10-7-10

Product Introduction

Our Midi series have been designed for routine horizontal gel electrophoresis with wide degrees of versatility. All of our gel tanks are made from injection molded construction for durable and leak proof environment. Varieties of special combs are available to maximize your sample size.

Midi plus-1 includes 10x7cm and 10x10cm gel tray to allow 100 samples per gel to be resolved. Ideal for those routinely check medium numbers of samples over short to medium gel run lengths.

FEATURES

- Low buffer consumption
- Easy sample loading
- Multi-channel pipette compatibility
- Great indications for gel making and running
- No tape, clamps or springs needed



Cat. No.	ME10-7-10
Dimension(mm)(WxLxH)	220x125x90
Gel dimension(mm)(WxL)	100x70 100x100
Maximum sample each tray(mm)	50 for 100x70 100 for 100x100
Rapid casting gel	Use casting dams

Cat. No.	ME10-7-10
Dimension (WxLxH)	8.7"x4.9"x3.5" (220x125x90mm)
Gel Dimension (WxL)	3.9"x2.8" (100x70mm) 3.9"x3.9" (100x100mm)
Maximum Sample	50 for 3.9"x2.8" (100x70mm) Tray (2x25 sample combs)
	100 for 3.9"x3.9" (100x100mm) Tray (4x25 sample combs)
Buffer Volume	0.08gal (300ml)
Construction	Injection molded construction
	Durable, leak proof environment for complete safety and long life
Cassette Type Electrode	Inexpensive, easy to replace
	Made of 99.99% corrosion resistant pure platinum
Electrical Safety	Lid can only be fastened in one way
	Upon lid removal, power is disconnected from buffer chamber
Rapid Casting Gel	Use casting dams