

## WAVER SHAKER

Waver shaker offers the optimal movement for the multiple staining and washing involved in gel, blotting and microscopy and immunostaining applications, its combination of vertical and horizontal orbital motion provides more rapid exchange of solvents and better mixing efficiency.

### Feature

- >Continuous or timed operation with automatic switch-off
- >Variable shaking speed from 5 to 100rpm for mixing and wash action
- >Two shaking moves in 2D or 3D
- >Adjustable platform angle
- >Interchangeable / stacking platforms, and accessories for a variety of vessels
- >15kg carry capability
- >One of 4 digital red LED display
- >High quality DC Brushless motor
- >Easy to extension shelf
- >Light weight for easy mobility
- >Very easy maintenance
- >33 x 33cm platform



## 2/3D Shaker

### Specification

### MW-23

Motion	Orbital and Nutation (2D or 3D)
Controller	Digital microprocessor controller
Speed / Inc	5 ~ 100rpm / 1rpm
Adjustable angle	0~25 °
Timer	1~9999 mins with alarm, continuous / 1 min
Carry Capability	15kg
Motor	DC Brushless type
Operation Temperature	Ambient to 40 °C
Platform Dimension (W x L)	33 x 33cm
Platform Material	Painted iron metal
Unit Dimension (W x L x H)	330 x 450 x 280
Construction	Painted iron metal
Weight	approx. 10.0kg
Rated Voltage	110V / 220V selectable

PROLAB



Proveedor de Laboratorios, S.A. de C.V.  
 Equipos y Reactivos para Investigación  
 Av. Siglo 21, # 50 Interior 116-A  
 Col. Santa Cruz de las Flores  
 Tlajomulco de Zúñiga, Jalisco 45640  
 Tel. 3338488484  
 email: [prolab@pro-lab.com.mx](mailto:prolab@pro-lab.com.mx)  
[www.prolab@pro-lab.com.mx](http://www.prolab@pro-lab.com.mx)