


PrO-Research Centrifuge K242 Large Prime



Features:

•PrO-Research range.

Designed specifically for Research, Universities, Pharmaceutical use.

W2t is standard to all and is a method normally found in Ultra Centrifuges

Some of the standard Quality features include: An extra thick stainless-steel bowl

- An alloy and steel frame
- A world leading inverter
- The best quality European brushless motor
- Our PrO-Research range of Centrifuges also include the following user features:
 - A Bright Blue LED display, 10 acceleration rates, 10 Deceleration rates, 108 memory Program. Program in Speed or Rcf (G) to 1 digit
 - A timer from 0-999 minutes, and Constant. Program in 1 second increments
 - A pulse short run button. W2t program
 - As centrifuges can be very dangerous in the event of a malfunction or if they are mishandled, we have included various safety features to protect the user. The Pro-Analytical centrifuge range includes the following safety features:
 - A multi-point locking lid
 - Emergency lid release
 - A lid spring strut
 - Lid lock detection
 - An over speed sensor
 - A barrier ring
 - A motor overheat sensor.



BRK5520 ROTOR
(4*500ml)

Model	BRK5520
Rotor type	4 x 500ml
Tube size	98 x 155 (85 x 130)mm
Speed mín, rpm	500
Speed máx, rpm	4000
Max. Acceleration, Rcf (G)	3400
Tube angle	0
Accel. time, seg	45 (35)
Desacel time, seg	45 (35)
Autoclave sterilization temp	121 (20) °C



Falcon with cap
Falcon with cap

point 15ml
point 50ml

17x120
29x115

ADAPTER	TUBES
AF815	48
AF850	20

Specifications:

Catalogue Nr.	PRO-K242
Speed	500-15,000 Rpm (1 Rpm steps)
RCF Max	22,000 G
Timer	0-9999 Min (1 sec steps)
Memory	108 programs
Acceleration rates	10 programs
Deceleration rates	10 programs
Rotor	BRK5520 * 1 (4*500ml)
Adapters	AF815*4, 12*15ml AF850*4, 5*50ml
Weight	63.5 Kg (without rotor)
Dimensions (H*W*D)	440*580*630mm
Power	750W
Voltage supply	120V, 50/60Hz

AF815 ADAPTER



AF850 ADAPTER

