



DONGSHENG BIOTECH

**PCR ENHANCER 500UL
(AUXILIAR DE AMPLIFICACION EN MUESTRAS DIFICILES)**



Description

PCR Enhancer System is used in conjunction with DNA polymerases to optimize PCR of problematic and/or GC-rich templates. The system offers higher primer specificity, broader magnesium concentration, broader annealing temperature, and improved thermostability DNA polymerases.

Usage

Add 2 – 8 µl PCR Enhancer to 25 µl amplification system, the other amplification conditions and procedure remain the same.

Store at -20°C